

# UrgoTul<sup>®</sup>

LIPIDO-COLLOID ULTRA FLEXIBLE CONTACT LAYER



CE 0459 

- Healing in a moist environment
- Painfree and atraumatic removal - does not dry out
- Proliferation of fibroblasts


## Indication:

Acute and chronic wounds, low to highly exuding, adapted to wounds located in awkward areas

## UrgoTul<sup>®</sup>

Supplied in boxes of individually pouched and sterile dressings, ready to use

Sizes	Units per box
5 x 5 cm	10
10 x 10 cm	10
15 x 20 cm	10
10 x 40 cm	10





# UrgoTul®



## DESCRIPTION

■ Urgotul is a sterile, non-adhesive, non occlusive lipido-colloid contact layer (issued from the Technology Lipido-Colloid, exclusive patented innovation by Laboratoires Urgo)

### COMPOSITION:

■ Urgotul is an ultra-flexible lipido-colloid dressing made of a polyester mesh impregnated with hydrocolloid (carboxymethylcellulose) and petroleum jelly particles

METHOD OF STERILIZATION: Sterilized by raditation.

## PROPERTIES

- On contact with wound exudate, hydrocolloid particles in Urgotul mass interact with petroleum jelly component to form a lipido-colloid interface creating a moist environment favourable for the healing process.
- Fatty in its chemical composition, without being greasy to the touch, Urgotul does not adhere to the wound or to the surrounding skin.
- Removal of the dressing is totally pain-free and atraumatic for the wound.
- Urgotul promotes fast healing by stimulating fibroblasts proliferation.
- Flexible and conformable, Urgotul is especially suitable for wounds located in awkward areas

## INDICATIONS

- Urgotul is indicated for the management of acute wounds (burns, skin tears, traumatic wounds, abrasions, post-operative wounds) and chronic wounds (leg ulcers, pressure ulcers and diabetic foot ulcers) at the granulation and epithelialisation stage.
- Urgotul is also intended for the management of epidermolysis bullosa.
- Ultra-flexible, Urgotul is particularly indicated to cover anfractuous or awkwardly placed wounds and to wick deep wounds.

## DIRECTIONS FOR USE

### METHOD OF USE

- Clean the wound as per local protocol.
- If an antiseptic is first used, carefully rinse the wound with normal sterile saline before applying UrgoTul.
- Dry the surrounding skin carefully.
- Remove the protective wings.
- Apply Urgotul to the wound.
- Urgotul can be cut with sterile equipment to fit the size of the wound.
- Cover Urgotul with a secondary dressing suitable to absorb wound exudate.
- Hold secondary dressing in place with a conforming bandage, an adhesive tape or an elasticated tubular bandage.
- Urgotul can be renewed every 2 to 4 days and left in place for up to 7 days depending on the wound condition.
- For patients with epidermolysis bullosa, Urgotul should be changed after 1 to 3 days.

### PRECAUTIONS

- Urgotul may adhere to latex gloves. In that case, moistening the gloves with saline solution will facilitate handling of the dressing.
- If clinical signs of local infection appear, the treatment may be continued under medical supervision, it may be preferable to use a Silver dressing instead.
- In case of deep wounds, Urgotul does not need to be applied in several layers.
- In the event of deep, anfractuous wounds or fistulas, leave part of the Urgotul dressing visible and accessible outside the wound.
- Do not re-sterilise the dressing.
- Store Urgotul dressings flat, protected from moisture and heat (temperature below 25°C).

# UrgoTul<sup>®</sup>

## Absorb Silicone Border

LIPIDO-COLLOID FOAM DRESSING  
WITH SILICONE ADHESIVE BORDER



CE 0459 

- Optimal healing result
- High absorption of exudate and drainage preventing any risk of maceration
- Waterproof silicone adhesive on the edges
- Pain-free and atraumatic care

### Indication:

Acute and chronic, moderately to heavily exuding wounds

### UrgoTul<sup>®</sup> Absorb Silicone Border

Supplied in boxes of individually pouched and sterile dressings, ready to use

Sizes	Units per box
8 x 8 cm	10
13 x 13 cm	10
15 x 20 cm	10
Sacrum: 20 x 20 cm	5





# UrgoTul<sup>®</sup>

## Absorb Silicone Border



### DESCRIPTION

- **UrgoTul Absorb Silicone Border** is an adhesive lipido-colloid foam dressing with soft-adherent TLC (Technology Lipido-Colloid), an innovative Technology developed by the Laboratoires Urgo.

### COMPOSITION:

- A soft-adherent TLC layer (polymers and hydrocolloid particles) combined with an absorbent polyurethane foam pad and a superabsorbent layer.
- A vapour permeable waterproof outer film with silicone adhesive on the edges.

**METHOD OF STERILIZATION:** Sterilized by ethylene oxide.

### PROPERTIES

- Soft-adherent TLC layer
- Optimal healing result
- Moist environment
- Pain-free dressing changes
- Fibroblasts proliferation
- Absorbent polyurethane foam pad and superabsorbent layer
- High fluid management, prevents maceration
- Suitable for use under compression bandaging
- Waterproof and silicone adhesive backing
- Soft and non-occlusive
- Very conformable to fit different wound shapes
- High vapour permeability to prevent maceration
- Good tolerance
- No need for secondary retention
- Mode of action
- In contact with wound exudate, TLC forms a gel and creates a moist wound environment, beneficial for the healing process enabling the key cells involved in the repair process (fibroblasts, keratinocytes and macrophages) to exert their action.

### INDICATIONS

- **UrgoTul Absorb Silicone Border** is indicated for moderately to heavily exuding acute wounds (such as burns, dermabrasions, traumatic wounds, postoperative wounds...) and chronic wounds (such as leg ulcers, pressure ulcers, diabetic foot ulcers).
- The sacrum version is recommended for exuding wounds located in the sacral area (sacral pressure ulcers...).

### DIRECTIONS FOR USE

#### METHOD OF USE

- Clean the wound as per local protocol.
- If an antiseptic is first used, carefully rinse the wound with normal sterile saline before applying **UrgoTul Absorb Silicone Border**.
- Dry the surrounding skin carefully.
- Remove the protective wings.
- Apply the soft-adherent side of **UrgoTul Absorb Silicone Border** to the wound (the adhesive silicone border should be at least 1 cm away from the wound).
- Smooth the dressing over the wound.
- Apply compression bandaging over the dressing when prescribed.
- To apply the sacrum version, place the dressing with the tip facing the sacral area.
- To remove **UrgoTul Absorb Silicone Border**, press down on the healthy skin, lift on the corner of the dressing and remove it carefully.
- **UrgoTul Absorb Silicone Border** dressing may be changed every 2 to 4 days, and left in place up to 7 days depending on the level of exudate and the condition of the wound.

#### PRECAUTIONS

- As it includes a super-absorbent layer, the central pad should not be cut. However, the adhesive edges can be cut if necessary using sterile scissors to fit the dressing to different wound shapes.
- If clinical signs of local infection are noted, treatment can be changed to an antibacterial dressing, dependent on clinical judgement.
- Any excess hair should be cut close to the skin to ensure good contact with the wound.
- In case of concomitant use with cream, an ointment, an emulsion, let the skin dry before the dressing application.
- Single use sterile individual packaging: re-using a single use dressing may lead to risks of infection.
- Do not re-sterilise the dressing.

#### CONTRAINDICATION

- Do not use if there is a known sensitivity to **UrgoTul Absorb Silicone Border** or its components.

# UrgoTul<sup>®</sup> Lite

LITE LIPIDO-COLLOID FOAM  
DRESSING



CE 0459 

- Healing in a moist environment
- Painfree and atraumatic removal
- Vertical absorption of exudate

## Indication:

Acute and chronic, low to moderately exuding wounds

## UrgoTul<sup>®</sup> Lite

Supplied in boxes of individually pouched and sterile dressings,  
ready to use

Sizes	Units per box
13 x 12 cm	16
15 x 20 cm	10





# UrgoTul<sup>®</sup> Lite



## DESCRIPTION

- UrgoTul Lite is a soft-adherent lipido-colloid foam dressing with TLC, an innovative Technology developed by the Laboratoires Urgo.

## COMPOSITION:

- A soft-adherent TLC layer made up of polymers and hydrocolloid particles (CarboxyMethylCellulose) combined with an absorbent polyurethane foam pad in contact with the wound,
- A protective non-woven polyurethane backing.

METHOD OF STERILIZATION: Sterilized by ionizing radiation

## PROPERTIES

- On contact with wound exudate, the hydrocolloid particles from the soft-adherent TLC layer form a lipido-colloid gel in contact with the wound.
- This gives the dressing specific properties:
  - creates a moist environment that promotes healing,
  - forms a gel in contact with the wound,
  - provides trauma-free removal without damaging any newly formed granulation tissue,
  - provides pain-free dressing changes for the patient,
  - respects the peri-wound skin.
- The absorbent polyurethane foam pad ensures drainage of exudate, protecting the peri-wound skin against maceration.
- The backing is soft, conformable and very comfortable: it helps ensure a good fit of the dressing to the wound shapes.
- The soft-adherent dressing facilitates nursing care regarding ease of use and patient compliance.

## INDICATIONS

- UrgoTul Lite is indicated for all low to moderately exuding chronic wounds (pressure ulcers, leg ulcers, diabetic foot ulcers), acute wounds (2nd degree burns, skin abrasions, traumatic wounds and post-operative wounds) and for the local treatment of fungating wounds in combination with a specific treatment.
- Due to the non-adhesive nature of the dressing, UrgoTul Lite is recommended in the treatment of wounds in which the surrounding skin is weakened.

## DIRECTIONS FOR USE

### METHOD OF USE

- Clean the wound as per local protocol, then rinse the wound with normal saline,
- Dry the surrounding skin carefully.
- UrgoTul Lite can be cut using sterile equipment to adjust the size of the dressing to fit the wound.
- Remove the protective wings,
- Apply the soft-adherent side of UrgoTul Lite to the wound,
- Secure the dressing in place with a suitable bandage,
- Apply compression bandaging over the dressing when prescribed.
- In the event of prior use of an antiseptic, carefully rinse the wound with normal sterile saline before applying UrgoTul Lite.
- UrgoTul Lite dressing may be changed 2 to 4 days, and left in place for up to 7 days depending on the exudate level and wound condition.

### PRECAUTIONS

- The soft-adherent TLC layer of UrgoTul Lite adheres to surgical gloves (latex). It is therefore recommended to avoid contact with the soft-adherent TLC layer.
- Due to the non-occlusive nature of this dressing, UrgoTul Lite can be used on infected wounds under medical supervision. However suitable treatment is recommended in colonised wounds prior to the use of UrgoTul Lite.
- Single use sterile individual packaging: re-using a single use dressing may lead to risks of infection.
- Do not re-sterilise.